




OPERATE AT 40 TO 60 KILOHERTZ.
RUN 10 TO 15 MINUTE CYCLES.
DEGAS SOLUTION BEFORE LOADING INSTRUMENTS.
KEEP THE UNIT TIGHTLY COVERED.

Ø
TOL

ACOUSTIC
CAVITATION ZONES
FREQUENCY RANGE:
40-60 kHz

WAVE
INTERACTION
POINT

CYCLE TIMER:
10-15 MIN RANGE



CAVITATION MECHANICALLY DISLODGES INVISIBLE PROTEIN RESIDUE FROM MICROSCOPIC CREVICES.

ULTRASONIC CLEANING PRECEDES STEAM STERILISATION.

TITLE:
ULTRASONIC CLEANING TANK CROSS-SECTION

SCALE: NTS
DATE: OCT 26, 2023
ENGINEER: J. DOE

DRAWING NO: UC-SCHEMATIC-001